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The screenshot displays the control software interface for a GC-MS system. The window title is "2000.40 - Not Ready". The interface is divided into several sections:

- Manual Control:** Includes buttons for "Auto Tune", "Shutdown", and "Reset".
- Control and Status:** Shows "Status: Idle" and "Vacuum System Pump Status: Ready".
- Temperatures:** Lists "Current Setpoints" for "Heated Zones":
  - Trap Temperature: 150 degrees C
  - Manifold Temperature: 75 degrees C
  - Transferline Temperature: 170 degrees CIt also shows "Vacuum System Pump: Full Speed".
- Operating Conditions:** Lists "Heated Zones":
  - Trap Temperature: 60 degrees C
  - Manifold Temperature: 74 degrees C
  - Transferline Temperature: 149 degrees CIt also shows "Vacuum System Turbo Speed: 99 %" and "Turbo Current: 2413 mA".

At the bottom, there is a status bar with the text: "Configuration : Aug 13 10:12:24 BOOTP Server is Waiting for a Request!". The Windows taskbar at the very bottom shows the time as 10:16 AM on August 13.

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The screenshot displays the 'System Control' software interface for a 2000.40 instrument. The interface is in Persian and shows various control panels and status indicators. The 'Diagnostics' tab is active, showing 'System Test' and 'Heater Test' options. The 'Monitor States' section shows 'Ion Gauge - On/Off' with 'Ion Gauge' checked and 'Filament 1' selected. The 'Vacuum System' shows 'Pump Status: Ready', 'Turbo Speed: 99%', and 'Turbo Current: 2414 mA'. The 'Heating System' shows 'Trap: 95', 'Manifold: 75', and 'Transferline: 174'. The 'Ion Gauge System' shows 'Vacuum Status: OK', 'Filament #1: OK', and 'Filament #2: OK'. The 'Reading' is 139.1 uTorr (Valid). The status bar at the bottom indicates 'Configuration : Aug 13 10:12:24 BOOTP Server is Waiting for a Request!'.

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The screenshot displays the 'System Control' software interface. The main window is titled '2000.40 - Not Ready'. The interface includes a menu bar (File, Edit, Inject, Automation, Recalculate, Instrument, Windows, Help) and a toolbar with various icons. The main control area is divided into several sections:

- Manual Control**: Includes 'Auto Tune', 'Start Bakeout', and 'Reset' buttons.
- Temperatures**: Contains 'Temperature (degrees C) Setpoints' and 'Bakeout Conditions' tables.
- Operating Conditions**: Displays current temperatures for Trap, Manifold, and Transferline.

Temperature (degrees C) Setpoints			
Trap	Manifold	Xferline	
150	75	170	

Bakeout Conditions			
Hold Time (h.)	Trap	Manifold	Xferline
20	210	110	170


Operating Conditions	
Trap Temperature:	85 degrees C
Manifold Temperature:	75 degrees C
Transferline Temperature:	171 degrees C


At the bottom of the window, a status bar reads: 'Configuration : Aug 13 10:12:24 BOOTP Server is Waiting for a Request!'. The Windows taskbar at the bottom shows the system clock as 10:19 AM on August 13.

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
Manual Control | Auto Tune | Temperatures | Diagnostics | Shutdown | Acquisition

Control and Status  
 State: Idle  
 Function: Start Auto Tune, Reset, Continue

Method: SetPoints  
 Air / Water Check  
 Electron Multiplier Tune  
 FC-43 Mass Calibration  
 Trap Function Calibration  
 Single Step

Operating Conditions

Hide Keypad | Spectrum and Event Message

---

Scan Number: 60, Time: 0.000 min. Base Peak: 330.1, Base Amount: 1  
 RIC: 237, Ion Range: 330.0 - 410.0 m/z Ion: 1000 us, Segment: 1, Channel: 1

10:23:00 Air/Water Check: Completed - Problems Found  
 10:23:00 Auto Tune: Completed  
 10:29:08 Auto Tune: Started  
 10:29:08 Integrator Zero Set: Started  
 10:29:22 Integrator Zero Set: Setting is OK (Setting: 30, Average Counts: 0.49)  
 10:29:22 Integrator Zero Set: Completed  
 10:29:22 Electron Multiplier: Started  
 10:29:24 Electron Multiplier: Offset Determined (Offset: 1 Count(s))  
 10:29:59 Electron Multiplier: Pre-Adjustment Successful (EM Voltage: 1700)

Configuration : Aug 13 10:12:24 BOOTP Server is Waiting for a Request!

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